



LARGE VOLUME BOTTLES

Multifunctional bottles made of PFA fluoropolymer

Ultra-
pure!



WHEN LARGE CAPACITY IS NEEDED

These bottles are an excellent choice when multifunctional use and superior material properties are required, which makes them perfect and suitable for a wide range of laboratory and industrial applications e.g. liquid handling in industrial processes, storage and shipping of both, aggressive media (acids) or sensitive active pharmaceutical ingredients (APIs).

FEATURES & BENEFITS

- ✓ Volumes from 3L up to 30L
- ✓ Contamination and loss-free handling of precious samples and APIs
- ✓ Usable for acids, pharmaceuticals and other media
- ✓ Multifunctional and flexible use
- ✓ Cleanable and reuseable
- ✓ FDA certification on request available

USAGE POSSIBILITIES



Media Supply



Media Collection



Storage



Transport

FEATURES



Translucent PFA

Rotomolded PFA



Various closures are available

- ✓ Made from ultra-pure PFA fluoropolymer
- ✓ Chemically inert
- ✓ Acid-resistant
- ✓ No interaction with bio-chemical substances
- ✓ Temperature-resistant (-200°C to +260°C)
- ✓ Unbreakable
- ✓ Autoclavable and reusable

APPLICATION SECTORS



Industrial processes



Pharmaceuticals



Packaging



Laboratory analysis

ORDERING INFORMATION

Order No.	Volume	Description	Neck ID	Material
T61-003	3 L	Narrow neck bottle with volume marks	37 mm	Translucent PFA
T61-005	5 L	Narrow neck bottle with volume marks	37 mm	Translucent PFA
T22-005	5 L	Narrow neck bottle with volume marks	33 mm	Rotomolded PFA (opaque)
T22-010	10 L	Narrow neck bottle with volume marks	33 mm	Rotomolded PFA (opaque)
T22-020	20 L	Narrow neck bottle with volume marks	33 mm	Rotomolded PFA (opaque)
T22-030	30 L	Narrow neck bottle with volume marks	33 mm	Rotomolded PFA (opaque)
T21-005	5 L	Wide neck bottle with volume marks	85 mm	Rotomolded PFA (opaque)
T21-010	10 L	Wide neck bottle with volume marks	99 mm	Rotomolded PFA (opaque)
T21-020	20 L	Wide neck bottle with volume marks	99 mm	Rotomolded PFA (opaque)
T21-030	30 L	Wide neck bottle with volume marks	99 mm	Rotomolded PFA (opaque)