

Guardian Autosampler Probe

Glass Expansion is pleased to release the Guardian autosampler probe. The unique design of the robust tip—which combines drip-resistance and built-in particle filtering—helps to prevent cross-contamination during probe movement and blockages in your nebulizer and capillary tubing.

Benefits

- Robust tip design prevents crushed and damaged tips due to misalignment.
- Drip-resistance prevents cross contamination of samples, especially with oils.
- Built-in particle filtering holds back particulates from blocking the line.
- Completely inert design, Ceramic, PEEK and PTFE construction.
- Interchangeable UniFit™ sample lines available in various IDs (e.g. 0.3, 0.50, 0.75 & 1.0mm)
- Designed to suit Teledyne Cetac®, Agilent®, PerkinElmer®, Shimadzu®, Aim Lab and Thermo Fisher Scientific™ Autosamplers.

Scan the below QR Code to watch the Guardian Autosampler Probe Performance Comparison video:



Please contact:
AHF analysentechnik AG
Kohlplattenweg 18
DE-72074 Tübingen, GERMANY
Tel.: +49 7071 53 952-00
info@ahf.de :: www.ahf.de

