


LedHUB® - High Power LED Light Engine



- Modular principle for up to 6 wavelengths between 365 and 840nm
- Field-upgradable
- Individual analogue and digital modulation inputs
- Precise intensity control in 0.1% steps (0-100%)
- Intelligent SYNC in- and outputs
- TEC temperature controlled LED chips
- Compatible to 2, 3 and 5mm Liquid Light Guides and SMA or FC/PC connectorized fibres
- Long lifetime – expected to exceed 25,000 hours of operating time



Available wavelengths:



Omicron-Laserage Laserprodukte GmbH

Lasers, Lightengines and LED lightsources for Life Science and Bio-Instrumentation



Distributor for LedHUB:

AHF analysentechnik AG
Kohlplattenweg 18
DE-72074 Tübingen, Germany

Tel.: +49 (0)7071 970 901-0
Fax: +49 (0)7071 970 901-99
info@ahf.de :: www.ahf.de