

Eluo Nebulizer Cleaning Tool

As easy as 1, 2, 3

1. Fill with solvent



2. Insert nebulizer



3. Clean nebulizer



Particle build-up in a nebulizer capillary and tip causes sample flow to be constricted, thereby reducing nebulizer efficiency and performance. Now, blocked nebulizers can be safely and easily restored to optimum performance with the Eluo.

The Eluo is designed to efficiently deliver a cleaning solution through the nebulizer capillary to dislodge particle build-up and thoroughly clean the nebulizer. One simple action does it all. No more messy procedures or shattered nebulizers in ultrasonic baths. Use the Eluo regularly to maintain nebulizer performance and prolong nebulizer life. Every lab should have an Eluo.

Order P/N [70-ELUO](#) for glass concentric nebulizers or P/N [70-ELUO-OPD](#) for OpalMist or DuraMist nebulizers.

Ordering Information

Part Number	Description
70-ELUO	Eluo Nebulizer Cleaning Tool (For glass concentric nebulizers)
70-ELUO-OPD	Eluo Nebulizer Cleaning Tool (For OpalMist or DuraMist nebulizers)



P/N 70-ELUO-OPD

For more information visit: www.geicp.com



GLASS EXPANSION
Quality By Design

Please contact:

AHF analysentechnik AG
Kohlplattenweg 18
72074 Tübingen
GERMANY



Tel.: +49 7071 53 952-00
Fax: +49 7071 53 952-99
info@ahf.de · www.ahf.de